



Physician Associates, PC

123 Main Street
Liverpool, MN 55555
(555) 123-4567

LOWER EXTREMITY ARTERIAL TESTING

Name Margaret Hughes	Telephone # (888) 123-4567	SSN/PIN	Sex Female	Age 79	Ht(in) 68	Date of Exam 05/01/2018
Referring Physician Don Holmberg	Telephone # (555) 123-4567	Account # 12344567890	Date of Birth 09/19/1938	Wt(lb) 170	Previous Exam 04/25/2017	

Risk Factors	Diabetes: No Hypertension: Yes Tobacco: Yes Hyperlipidemia: Yes Body Mass Index: 25.8 (overweight)
Assoc. Illness	Cardiac: Yes Carotid: Unknown Pulmonary: Unknown Renal: Unknown
Indications	Claudication (both sides)
Notes	Patient is currently on medication for Hypertension.

Right	Results	Left
<p>Pressures</p> <p>Brachial 114</p> <p>Above Knee 119 SBI 1.04</p> <p>Calf 97 SBI 0.85</p> <p>Ankle 74 ABI 0.65</p> <p>Toe 17 TBI 0.15</p> <p>PVR Waveforms</p> <p>Above Knee Normal</p> <p>Calf Mild</p> <p>Ankle Moderate</p> <p>Toe Severe</p>	<p>ABI RANGE VALUE*</p> <p>>1.30 Noncompressible</p> <p>1.00 - 1.29 Normal</p> <p>0.91 - 0.99 Borderline (equivocal)</p> <p>0.41 - 0.90 Mild to moderate</p> <p>0.00 - 0.40 Severe P.A.D.</p> <p>*P.A.D. Coalition *N Engl J Med 2001; 344: 1608-21</p>	<p>Pressures</p> <p>Brachial 114</p> <p>Above Knee 119 SBI 1.04</p> <p>Calf 97 SBI 0.85</p> <p>Ankle 74 ABI 0.65</p> <p>Toe 17 TBI 0.15</p> <p>PVR Waveforms</p> <p>Above Knee Normal</p> <p>Calf Mild</p> <p>Ankle Moderate</p> <p>Toe Severe</p>

Interpretation

This is the area where the professional interpretation is included.

Additional Comments

This is where additional notes can be added.

Dee Brown

Electronic signature

Dee Brown
Technologist: 05/02/2018


William Howard

Electronic signature

William Howard
Provider: 05/03/2018

All reports are for illustrative purposes only. Actual appearance may vary.





Ankle systolic blood pressure

Pulse Volume Recording (PVR) waveforms at above knee, calf, ankle and toe

Arm systolic blood pressure

Gain setting affects height (amplitude) of waveforms

Toe systolic blood pressure

SBI values

ABI value

TBI value

All reports are for illustrative purposes only. Actual appearance may vary.

